

Shortest CEP stabilised Ti:Sapphire pulse

Abstract: A dual-output, self-synchronised, CEP-stabilised Ti:sapphire oscillator is presented, with a pulse duration of 4.4fs and a CEP phase noise of < 100 mrad using a Finesse 6 Pure pump laser.

Since the first demonstration of the gain material Ti:sapphire as a tuneable laser medium, there has been much interest in the exploitation of its broadband optical output, chief among which is the generation of ultra-short optical pulses. Since the pulse duration is inversely proportional to the optical bandwidth supported by the Ti:sapphire oscillator, careful attention must be paid to the optical design to compensate for dispersive effects.

The generation of sub-20 femtoseconds pulses (fs or 10^{-15} seconds) is now routinely demonstrated and has been commercially available for several years. In recent years, however, the pulse duration has been driven down to sub-10 fs by exploiting the full bandwidth of the Ti:sapphire crystal. Once the pulse duration is reduced to sub-6 fs, there are only approximately two cycles of the electric field within the envelope of the optical pulse. Figure 1 shows a comparison of two femtosecond pulses, their pulse envelope plotted with the corresponding electric field oscillation.




Figure 1: A comparison of two femtosecond pulses with the corresponding electric field oscillations. In the few-cycles pulse shown, the Ti:sapphire centre frequency of 375 THz, a wavelength of 800 nm respectively, the period of the electric field oscillation is 2.7 fs.




Figure 2: An image of the Finesse 6 Pure (Laser Quantum Ltd) pumping an octave-spanning VENTEON | PULSE : ONE femtosecond oscillator

It can be seen that the number of contributing field oscillations supported by the longer pulse - which is already as short as 15 fs - is significantly higher than for the few-cycle femtosecond pulse. By stabilizing the carrier envelope phase (CEP, ϕ_0), the electric field of such few-cycle pulses can be controlled opening up for novel, significant ultrafast experiments such as high harmonic generation and attosecond science.

Using a Finesse Pure 6 W pump source at 532 nm from Laser Quantum Ltd, integrated into a common water-cooled platform, VENTEON Laser Technologies GmbH has demonstrated an octave-spanning direct CEP-stabilised output of well below 5 fs without the need for additional spectral broadening (Figure 2).

This system delivers the shortest commercially available pulses directly from a Ti:sapphire oscillator, with a SPIDER characterised pulse as short as 4.4 fs as shown in Figure 3.

The CEP phase noise of those pulses is less than 100 mrad. This requires a minimal amplitude noise of the Finesse Pure within the whole noise spectrum. With RMS noise significantly less than 0.03% (10 Hz to 100 MHz), the Finesse Pure is ideally suited to this purpose.

This laser system is directly applicable for optical parametric chirped pulse amplification (OPCPA) applications, since its octave-spanning

spectral bandwidth allows the system to be equipped with two separate outputs which are intrinsically self-synchronised with ultra-low timing jitter. The first output can provide CEP-stable broadband sub-6 fs signal pulses; the second narrowband output at a centre frequency of 1030 nm can deliver a pulse energy of > 30 pJ in a spectral bandwidth of 10 nm (FWHM) and is an ideal seed for Yb-doped amplifier stages used for pumping NOPA stages of an OPCPA amplifier system.

In such an OPCPA amplifier system, the ultrashort pulse duration of the sub-6 fs signal seed pulses can be preserved and amplified up to multi- μ J pulse energies making an ideal source for high harmonic generation and ultrafast spectroscopy.




Figure 3: A SPIDER characterised pulse of 4.4 fs is shown: the shortest commercially available pulse from a CEP-stabilised Ti:sapphire laser.

Acknowledgements: Dr Stefan Rausch and Dr Thomas Binhammer of VENTEON Laser Technologies GmbH and Dr Lawrie Gloster of Laser Quantum Ltd.